

**AUR-200 URINE ANALYZER**

AUR-200 Urine Analyzer is a device designed for early screening and detecting kidney diseases, diabetes, urinary tract infection and liver diseases.

**Advantages:**

- ◆ Comprehensive Range of Test Menu
- ◆ Built-in Battery Available for STAT Usage
- ◆ Testing Up to 14 Parameters in One Go
- ◆ Light Weight for Easy Mobility
- ◆ Built-in Thermal Printer for Instant Printing
- ◆ Touch Screen Operation with Big Screen for Easy Operation

**TECHNICAL SPECIFICATIONS**

Model	AUR-200
Principle	Photoelectric Colorimetry
Testing Speed	60 tests/hour
Data Memory	10.000 results
Strip Incubation Time	1 minute
Light Wavelength	470nm, 525nm and 625nm
Display	3.5" LCD touch screen
Printer	Built-in thermal printer
Connectivity	RS232 port, USB 2 port
Working Condition	10-30°C, RH≤80%, 76-106kPa
Working Voltage	AC100-240V, 50/60Hz, DC 5.85V±0.45V (AAx4)
Weight	600g
Dimension	200*137*52mm