


## TDZ4-WS TABLETOP LOW SPEED CENTRIFUGE

TDZ4-WS Tabletop Low Speed Centrifuge with its capacity of 18x10ml, 12x20ml and 4x50ml, which it is widely used in most of the hospitals and universities by its high performance.

### Features:

- ◆ Microprocessor control, DC brushless motor.
- ◆ Touch panel and LCD display.
- ◆ Special damping structure to reduce vibration.
- ◆ Whole steel housing.
- ◆ Countdown time showed by seconds if the centrifugation time is less than 1min.



### TECHNICAL PARAMETER

Model	TDZ4-WS
Max. Speed	4000rpm
Max. RCF	2200xg
Max. Capacity	12 x 20ml
Speed Accuracy	±30rpm
Time Setting Range	1min to 99min
Noise	<62dB(A)
Power Supply	AC220V±22V 50/60Hz 5A
Total Power	100W
Dimension (W x D x H)	330x390x250mm
Package Dimension (W x D x H)	395x475x290mm
Net Weight	16kg



Capacity: 12 x 20ml  
Max. Speed: 4000rpm  
Max. RCF: 2200xg



Capacity: 18 x 10ml  
Max. Speed: 4000rpm  
Max. RCF: 2200xg



Capacity: 4x50ml  
Max. Speed: 4000rpm  
Max. RCF: 2200xg