

H1650 TABLETOP HIGH SPEED CENTRIFUGE

- ◆ Microprocessor control, DC brushless frequency conversion motor.
- ◆ Operated by touch panel, LCD display.
- ◆ The electronic door interlock for safe operation and stainless steel chamber.
- ◆ It is programmable and can store at least ten programs in memory.
- ◆ RCF value can be calculated automatically.
- ◆ Many kinds of rotors are available for convenient use.



TECHNICAL PARAMETER

Model	H1650
Max. Speed	16500rpm
Max. RCF	18930xg
Max. Capacity	12x10ml
Speed Accuracy	±50rpm
Time Setting Range	1min to 99min59s
Noise	≤65dB(A)
Power Supply	AC220V±22V 50/60Hz 5A
Total Power	400w
Dimension (W x D x H)	325x420x325mm
Package Dimension(W x D x H)	385x640x330mm
Net Weight	20kg



Capacity: 12x1.5/2.2ml
Max. Speed: 16500rpm
Max. RCF: 18930xg



Capacity: 12x5ml
Max. Speed: 13500rpm
Max. RCF: 12835xg



Capacity: 24x1.5/2.2ml
Max. Speed: 13500rpm
Max. RCF: 16280xg



Capacity: 48x0.5ml
Max. Speed: 13500rpm
Max. RCF: 13810xg



Capacity: 12x10ml
Max. Speed: 12300rpm
Max. RCF: 15550xg



Capacity: 4x10ml
Max. Speed: 6000rpm
Max. RCF: 3200xg



Capacity: 2x2x48well
Max. Speed: 3000rpm
Max. RCF: 1400xg

